



EMAL6 Rabbit pAb

货号: **AYT1115**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB 1:500-2000:
理论分子量	215kD
实测分子量	215kD
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Cytoplasm , cytoskeleton .
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from human EMAL6
------	--

靶点信息

研究背景	The cytoskeleton is a dynamic three-dimensional filamentous structure in the cytoplasm. Its roles include maintenance of cell shape, cell movement (in eukaryotes), cytokinesis, and the organization of organelles or organelle-like structures within the cell. The cytoskeleton includes microfilaments (actin-like proteins), microtubules (tubulin-like proteins), the intermediate filaments (mostly in eukaryotes) and the MinD-ParA proteins, which appear to be unique to prokaryotes.
基因ID	400954;
基因名	EML6 EML5L
Swiss	Q6ZMW3 (https://www.uniprot.org/uniprotkb/Q6ZMW3/entry)
别名	EML6 EML5L, Echinoderm microtubule-associated protein-like 6 (EMAP-6) (Echinoderm microtubule-associated protein-like 5-like)

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)