



DAPK3 Rabbit mAb

货号: **AYT1121**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC ICC/IF
推荐浓度	WB: 1:1000 - 1:5000 IHC: 1:100- 1:300 IF/ICC: 1:100 - 1:600
理论分子量	53kDa
实测分子量	54kD
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	mouse brain,mouse lung,rat brain
细胞定位	Cytoplasm,Cytoplasm,Nucleus
纯化	Affinity purification


抗原信息

抗原信息	Synthesized peptide derived from Human ZIP Kinase.
------	--

靶点信息

研究背景	Death-associated protein kinase 3 (DAPK3) induces morphological changes in apoptosis when overexpressed in mammalian cells. These results suggest that DAPK3 may play a role in the induction of apoptosis.
基因ID	1613
基因名	DAPK3
Swiss	O43293 (https://www.uniprot.org/uniprotkb/O43293/entry)
别名	DAPK3,DLK,ZIP,ZIPK,ZIP Kinase Antibody,DAP-like kinase,MYPT1 kinase,Zipper-interacting protein kinase, ZIP Kinase

产品验证



Immunofluorescence analysis of A431 cells labelling DAPK3 with 1:200 . The cells were fixed with cold 10 0% methanol (10min, 4°C) and blocked in 1% BSA/10% normal goat serum/0.3M glycine in 0.1% PBS-Twe en 20 for 1h. The cells were then incubate with Goat Anti-Rabbit IgG Fluor488 (1:200, shown in green) at room temprature for 1h.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)