



PJA2 Rabbit pAb

货号: **AYT1123**

产品信息

反应	
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC ICC/IF
推荐浓度	WB 1:500-2000;: IHC 1:50-300: IF 1:50-200:
理论分子量	85kD
实测分子量	100-180 kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.324
偶联物	Unconjugated
阳性对照	
细胞定位	Cytoplasm . Cell membrane
纯化	亲和纯化

抗原信息



抗原信息	The antiserum was produced against synthesized peptide derived from human PJA2.
------	---

靶点信息


研究背景	function:Has E2-dependent E3 ubiquitin-protein ligase activity. ,pathway:Protein modification; protein ubiquitination. ,similarity:Contains 1 RING-type zinc finger. ,subcellular location:Localizes at the cytoplasmic side of endoplasmic reticulum and Golgi apparatus. Expressed in the postsynaptic density region of synapses. ,subunit:Binds ubiquitin-conjugating enzymes (E2s) . In vitro , interacts with the ubiquitin-conjugating enzyme , UBE2D2.
基因ID	9867;
基因名	PJA2
Swiss	O43164 (https://www.uniprot.org/uniprotkb/O43164/entry)
别名	KIAA0438; Neurodap 1; Neurodap1; PJA 2; PJA2; PJA2_HUMAN; Praja 2; Praja 2 RING H2 motif containing; praja ring finger 2; Praja ring finger 2 E3 ubiquitin protein ligase; Praja2; RING finger protein 131; RNF 131; RNF131;

产品验证

Western blot analysis of extracts from K562 cells using PJA2 at 1:500.



Human cerebral cortex was stained with anti-PJA2 rabbit antibody, 1: 100



Human Liver was stained with anti-PJA2 rabbit antibody, 1: 100

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)