



ARHGAP32 Rabbit pAb

货号: **AYT1132**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB : 1:500-1:2000
理论分子量	220 kDa
实测分子量	231 kDa
形式	Liquid
保存条件	Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	
阳性对照	
细胞定位	
纯化	亲和纯化

抗原信息

抗原信息	Synthesized peptide derived from human ARHGAP32
------	---


靶点信息

研究背景	ARHGAP32, also known as RICS, belongs to the PX domain-containing GAP family. ARHGAP32 encodes a neuron-associated GTPase-activating protein that regulates dendritic spine morphology and strength. ARHGAP32 also can interact with RhoA and β -catenin, thereby inhibiting phosphorylation of LIMK/cofilin and promoting epithelial-to-mesenchymal transition (EMT) process. ARHGAP32 is highly expressed in brain and testis
------	---

基因ID	9743
基因名	ARHGAP32
Swiss	A7KAX9 (https://www.uniprot.org/uniprotkb/A7KAX9/entry)
别名	ARHGAP32, GC GAP, GRIT, KIAA0712, p200RhoGAP, p250GAP, Rho GTPase activating protein, RICS

产品验证

Western blot analysis of extracts from A431 cell using ARHGAP32 at 1:500.



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)