



GAP43 Mouse mAb

货号: **AYT1135**

产品信息

反应	Human,Mouse,Rat
宿主	Mouse
克隆性	Monoclonal
预测反应	
应用	WB IHC ICC/IF
推荐浓度	WB: 1:500-2000 I HC : HC 1:100-500; IF 1:50-200:
理论分子量	25kD
实测分子量	35-43kD
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Cell membrane>Peripheral membrane protein
纯化	亲和纯化


抗原信息

抗原信息	Recombinant Protein of GAP-43
------	-------------------------------


靶点信息

研究背景	The protein encoded by this gene has been termed a 'growth' or 'plasticity' protein because it is expressed at high levels in neuronal growth cones during development and axonal regeneration. This protein is considered a crucial component of an effective regenerative response in the nervous system. Alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. [provided by RefSeq
基因ID	2596;
基因名	GAP43
Swiss	P17677 (https://www.uniprot.org/uniprotkb/P17677/entry)
别名	GAP43 ; Neuromodulin ; Axonal membrane protein GAP-43 ;

产品验证



Western blot analysis of extracts from 3T3 cells and Mouse brain tissue using GAP43 at 1:500.



Immunohistochemical analysis of paraffin-embedded Rat Brain Tissue using GAP-43 Monoclonal Antibody., 1:100

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)