



G-CSFR Rabbit pAb

货号: **AYT1137**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB 1:500-1:2000: IHC : 1:50-1:200
理论分子量	95kD
实测分子量	70-95kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Cell membrane ; Single-pass type I membrane protein
纯化	亲和纯化

抗原信息

抗原信息	The antiserum was produced against synthesized peptide derived from the Internal region of human CSF 3R
------	---

靶点信息

研究背景	The protein encoded by this gene is the receptor for colony stimulating factor 3 , a cytokine that controls the production , differentiation , and function of granulocytes. The encoded protein , which is a member o f the family of cytokine receptors , may also function in some cell surface adhesion or recognition process es. Alternatively spliced transcript variants have been described. Mutations in this gene are a cause of Ko stmann syndrome , also known as severe congenital neutropenia
基因ID	1441;
基因名	CSF3R
Swiss	Q99062 (https://www.uniprot.org/uniprotkb/Q99062/entry)
别名	CSF3R ; GCSFR ; Granulocyte colony-stimulating factor receptor ; G-CSF receptor ; G-CSF-R ; CD114

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)