



DSCAM Rabbit pAb

货号: **AYT1144**

产品信息

反应	
宿主	
克隆性	
预测反应	
应用	
推荐浓度	
理论分子量	222kD
实测分子量	230KD-270 kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Secreted. Cell membrane>Single-pass type I membrane protein. Cell projection>Axon. Cell projection>Dendrite. Cell projection>Growth cone. Cell junction>Synapse.
纯化	Affinity purification

抗原信息


抗原信息	Synthesized peptide derived from part region of human protein.
------	--

靶点信息

研究背景	This gene is a member of the immunoglobulin superfamily of cell adhesion molecules (Ig-CAMs) , and is involved in human central and peripheral nervous system development. This gene is a candidate for Down syndrome and congenital heart disease (DSCHD) . A gene encoding a similar Ig-CAM protein is located on chromosome 11. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene
基因ID	1826;
基因名	DSCAM
Swiss	O60469 (https://www.uniprot.org/uniprotkb/O60469/entry)
别名	CHD2 42; CHD2 52; CHD2; Down syndrome cell adhesion molecule; DSCAM; DSCAM_HUMAN;

产品验证

Western blot analysis of extracts from C6 cells using DSCAM at 1:500.



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)