

ARID2 Rabbit pAb

货号: **B10294**

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens , Mus musculus
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	191kDa/197kDa
实测分子量	250kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	LO2
细胞定位	Nucleus
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1586-1835 of human ARID2 (NP_689854.2).
序列	QNVVPQNTPMPPSPAVQVQGQPNSSQSPFSGSSQPGDPMRKPGQNFMCCLWQSCKKWFQTPSQVFYHAATEHGGKD VYPGQCLWEGCEPFQQRFSFITHLQDKHCKSDALLAGLKQDEPGQAGSQKSSTKQPTVGGTSSTPRAQKAIVNHPSAA LMALRRGSRNLVFRDFTDEKEGPITKHIRLTAALILKNIGKYSECGRLLKRHENNLVLAISNMEASSTLAKCLYELNFTVQS KEQEKDSEMLQ

靶点信息

研究背景	This gene encodes a member of the AT-rich interactive domain (ARID)-containing family of DNA-binding proteins. Members of the ARID family have roles in embryonic patterning, cell lineage gene regulation, cell cycle control, transcriptional regulation and chromatin structure modification. This protein functions as a subunit of the polybromo- and BRG1-associated factor or PBAF (SWI/SNF-B) chromatin remodeling complex which facilitates ligand-dependent transcriptional activation by nuclear receptors. Mutations in this gene are associated with hepatocellular carcinomas. A pseudogene of this gene is found on chromosome1. [provided by RefSeq, Dec 2016]
基因ID	196528
基因名	ARID2
Swiss	Q68CP9
别名	ARID2;BAF200;p200

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn