

IFITM3 Rabbit pAb

货号: **B10316**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB IHC: Homo sapiens WB: Homo sapiens
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:100
理论分子量	14kDa
实测分子量	15kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa
细胞定位	Cell membrane,Late endosome membrane,Lysosome membrane,Single-pass type II membrane protein
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-133 of human IFITM3 (NP_066362.2).
序列	MNHTVQTFFSPVNSGQPPNYEMLKEEHEVAVLGAPHNPAPTSTVIHIRSETSVDPDHVVWSLFNTLFMNPCCLGFI AFAYS VKSRDRKMGVDVTGAQAYASTAKCLNIWALILGILMTILLIVIPVLIFQAYG

靶点信息

研究背景	The protein encoded by this gene is an interferon-induced membrane protein that helps confer immunity to influenza A H1N1 virus, West Nile virus, and dengue virus. Two transcript variants, only one of them pr otein-coding, have been found for this gene. Another variant encoding an N-terminally truncated isoform has been reported, but the full-length nature of this variant has not been determined.
------	---

基因ID	10410
基因名	IFITM3
Swiss	Q01628
别名	IFITM3;1-8U;DSPA2b;IP15

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn