

RAB5A Rabbit pAb

货号: **B10317**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	22kDa/23kDa
实测分子量	24KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,Mouse brain,Rat brain
细胞定位	Cell membrane,Cell projection,Cytoplasm,Cytoplasmic side,Cytoplasmic vesicle,Early endosome membrane,Lipid-anchor,Melanosome,Membrane,cytosol,ruffle
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-215 of human RAB5A (NP_004153.2).
序列	MASRGATRPNGPNTGNKICQFKLVLLGESAVGKSSLVLRFBVKGQFHEFQESTIGAAFLTQTVCLDDTTVKFEIWDTAGQERYHSLAPMYRGAQAIVVYDITNEESFARAKNWVKELQRQASPNIVIALSGNKADLANKRAVDQEAQSYADDNSLLFMETSAKTSMNVNEIFMAIAKKLPKNEPQNPGANSARGRGVDLTEPTQPTRNQCCSN

靶点信息

研究背景	RAB5A (RAB5A, member RAS oncogene family) is a protein-coding gene. Diseases associated with RAB5A include borna disease, and venezuelan equine encephalitis. GO annotations related to this gene include G DP binding and GTP binding. An important paralog of this gene is RAB20.
------	---

基因ID	5868
基因名	RAB5A
Swiss	P20339
别名	RAB5A;RAB5

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn