

ANXA1 Rabbit pAb

货号: **B10323**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Huh7.5.1,HepG2 , Homo sapiens
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	38kDa
实测分子量	38KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	C2C12,C6
细胞定位	Apical cell membrane,Basolateral cell membrane,Cell membrane,Cell projection,Cytoplasm,Cytoplasmic vesicle,Cytoplasmic vesicle membrane,Early endosome,Endosome membrane,Extracellular side,Lateral cell membrane,Membrane,Nucleus,Peripheral membrane protein,Secreted,cilium,exosome,extracellular space,phagocytic cup,secretory vesicle lumen
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-346 of human ANXA1 (NP_000691.1).
序列	MAMVSEFLKQAWFIENEEQEYVQTVKSSKGGPGSAVSPYPTFNPSSDVAALHKAIMVKGVDIATIDILTKRNNAQRQQIK AAYLQETGKPLDETLKKALTGHLEEVVLALLKTPAQFDADELRAAMKGLGTDEDTLIEILASRTNKEIRDINRVYREELKRD KDITSDTSGDFRNALLSLAKGDRSEDFGVNEDLADSDARALYEAGERRKGTDVNVFNTILTTRSYPQLRRVFQKYTKYSK HDMNKVLDLELKGDIKCLTAIVKCATSKPAFFAEKLHQAMKGVGTRHKALIRIMVSRSEIDMNDIKAFYQKMYGISLCQAIL DETKGDYEKILVALCGGN

靶点信息

研究背景	This gene encodes a membrane-localized protein that binds phospholipids. This protein inhibits phospholipase A2 and has anti-inflammatory activity. Loss of function or expression of this gene has been detected in multiple tumors.
基因ID	301
基因名	ANXA1
Swiss	P04083
别名	ANXA1;ANX1;LPC1

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn