

ARHGEF1 Rabbit pAb

货号: **B10340**

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:100 IF/ICC: 1:50 - 1:200
理论分子量	98kDa/102kDa/104kDa/105kDa
实测分子量	115kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Jurkat,LO2
细胞定位	Cytoplasm,Membrane
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 783-912 of human ARH GEF1 (NP_004697.2).
序列	PLLSSSENGNGGRETS PADARTERILSDLLPFCRPGPEGQLAATALRKVLSLKQLLFP AEEDNGAGPPRDGDGVPGGGPL SPARTQEIQENLLSLEETMKQLEEL EEEFCRLRPLLSQLGGNSVPQPGCT

靶点信息

研究背景	Rho GTPases play a fundamental role in numerous cellular processes that are initiated by extracellular stimuli that work through G protein coupled receptors. The encoded protein may form complex with G proteins and stimulate Rho-dependent signals. Multiple alternatively spliced transcript variants have been found for this gene, but the full-length nature of some variants has not been defined.
------	---

基因ID	9138
基因名	ARHGEF1
Swiss	Q92888
别名	ARHGEF1;GEF1;LBCL2;LSC;P115-RHOGEF;SUB1.5

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn