

PFKP Rabbit pAb

货号: **B10346**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:100 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	85kDa
实测分子量	86KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,Jurkat,Mouse brain,Mouse heart,Rat testis,Rat brain
细胞定位	Cytoplasm
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 336-540 of human PFKP (NP_002618.1).
序列	ATPDTPACVVSLNGNHAVRLPLMECVQMTQDVQKAMDERRFQDAVRLRGRSFAGNLNTYKRLAIKLPDDQIPKTNCNVA VINVGAPAAGMNAAVRSAVRVGIADGHRMLAIYDGFDFGAKGQIKEIGWTDVGGWTGQGGSGILGTRVLPKYLEEIAT QMRTHSINALLIIGGFEAYLGLELSAAREKHFEFCVPMVMVPATVS

靶点信息

研究背景	This gene encodes a member of the phosphofructokinase A protein family. The encoded enzyme is the pl atelet-specific isoform of phosphofructokinase and plays a key role in glycolysis regulation. This gene ma y play a role in metabolic reprogramming in some cancers, including clear cell renal cell carcinomas, and cancer of the bladder, breast, and lung. Alternative splicing results in multiple transcript variants.
------	---

基因ID	5214
基因名	PFKP
Swiss	Q01813
别名	PFKP;ATP-PFK;PFK-C;PFK-P;PFKF

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn