

CAPG Rabbit pAb

货号: **B10349**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	36kDa/38kDa
实测分子量	38kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	THP-1,HL-60,Mouse lung,Mouse spleen,Rat heart,Rat lung
细胞定位	Cytoplasm,Melanosome,Nucleus
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-348 of human CAPG (NP_001738.2).
序列	MYTAIPQSGSPFPGSVQDPGLHVWRVEKLKVPVQAQENQGVFFSGDSYLVLHNGPEEVSHLHLWIGQQSSRDEQGACA VLAVHLNTLLGERPVQHREVQGNESDLFMSYFPRGLKYQEGGVESAFHKTSTGAPAAIKKLYQVKGKKNIRATERALNWD SFNTGDCFILDGQNIFAWCGGKSNILERNKARDLALAIRDSEKQGAQVEIVTDGEEPAEMIQVLGPKPALKEGNPEEDLT ADKANAQAAAALYKVSDATGQMNLTQVADSSPFALELLISDDCFVLDNGLCGKIYWKGRKANEKERQAALQVAEGFISRM QYAPNTQVEILPQGHESPIFKQFFKDWK

靶点信息

研究背景	This gene encodes a member of the gelsolin/villin family of actin-regulatory proteins. The encoded protein reversibly blocks the barbed ends of F-actin filaments in a Ca ²⁺ and phosphoinositide-regulated manner, but does not sever preformed actin filaments. By capping the barbed ends of actin filaments, the encoded protein contributes to the control of actin-based motility in non-muscle cells. Alternatively spliced transcript variants have been observed for this gene.
基因ID	822
基因名	CAPG
Swiss	P40121
别名	CAPG;AFCP;HEL-S-66;MCP

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn