

TAGLN2 Rabbit pAb

货号: **B10350**

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	22kDa/24kDa
实测分子量	22KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,A-549,Mouse heart
细胞定位	cytoskeleton,cytosol,extracellular exosome,extracellular region
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-199 of human TAGLN 2 (NP_003555.1).
序列	MANRGPAYGLSREVQQKIEKQYDADLEQILIQWITTQCRKDVGRPQPGRENFQNLKDGTVLCELINALYPEGQAPVKKI QASTMAFKQMEQISQFLQAAERYGINTTDIFQTVDLWEGKNMACVQRTLMLNLGGLAVARDDGLFSGDPNWFPPKSKEN PRNFSDNQLQEGKNVIGLQMGNTNRGASQAGMTGYGMPRQIL

靶点信息

研究背景	The protein encoded by this gene is similar to the protein transgelin, which is one of the earliest markers of differentiated smooth muscle. The specific function of this protein has not yet been determined, although it is thought to be a tumor suppressor. Multiple transcript variants encoding different isoforms have been found for this gene.
------	--

基因ID	8407
基因名	TAGLN2
Swiss	P37802
别名	TAGLN2;HA1756

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn