

# CRKL Rabbit pAb

货号: B10367

## 产品信息

|       |                                                                                                      |
|-------|------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                                      |
| 宿主    | Rabbit                                                                                               |
| 克隆性   | Polyclonal                                                                                           |
| 预测反应  |                                                                                                      |
| 应用    | WB IF/ICC                                                                                            |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000<br><b>IF/ICC:</b> 1:50 - 1:200                                             |
| 理论分子量 | 33kDa                                                                                                |
| 实测分子量 | 37kDa                                                                                                |
| 形式    | Liquid                                                                                               |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                         |
| 阳性对照  | LO2,Jurkat,SKOV3,Mouse brain,Rat heart                                                               |
| 细胞定位  | cytosol,extrinsic component of postsynaptic membrane,neuromuscular junction,nucleoplasm              |
| 纯化    | Affinity purification                                                                                |

## 抗原信息

|      |                                                                                                                                                                                                                                                                                                                              |
|------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 1-303 of human CRKL ( NP_005198.1).                                                                                                                                                                                                            |
| 序列   | MSSARFDSSDRSAWYMGPVSRQEAQTRLQGQRHGMFLVRDSSTCPGDYVLSVSENSRVSHYIINSLPNRRFKIGDQEFD<br>HLPALLEFYKIHLYLDTTTLIEPAPRYSPPPMGSVSAPNLPTAEDNLEYVRTLYDFPGNDAEDLPFKKGEILVIIKPEEQWWSA<br>RNKDGRVGMIPVPYVEKLVRSPPHGKHGHRNSNSYGIPEPAHAYAQPQTTTLPVAVSGSPGAAITPLPSTQNGPVPFAKAIQ<br>KRVPCAYDKTALALEVGDIVKVTRMNINGQWEGEVNGRKGLFPFTHVKIFDPQNPDENE |

## 靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                 |
|------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This gene encodes a protein kinase containing SH2 and SH3 (src homology) domains which has been shown to activate the RAS and JUN kinase signaling pathways and transform fibroblasts in a RAS-dependent fashion. It is a substrate of the BCR-ABL tyrosine kinase, plays a role in fibroblast transformation by BCR-ABL, and may be oncogenic. |
|------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|       |        |
|-------|--------|
| 基因ID  | 1399   |
| 基因名   | CRKL   |
| Swiss | P46109 |
| 别名    | CRKL   |

产品验证

实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)