

PARG Rabbit pAb

货号: B10436

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:100 - 1:500
理论分子量	54kDa/60kDa/99kDa/102kDa/111kDa
实测分子量	111KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	PC-3
细胞定位	Cytoplasm,Mitochondrion,Mitochondrion matrix,Nucleus
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 777-976 of human PARG (NP_003622.2).
序列	HNECLIITGTEQYSEYTGyaETyrWSRSHEDGserDDWQRRCTEIVaIDALHFRRyLDQFVPEKMRRELNKAYCGFLRPgV SSEnLSAVATGNWGCGAFGGDARLKALiQILAAAAAERDVVYFTFGDSELMRDIYSMHIFLTERKLTVGdVYKLLLRyYNEE CRNCSTPGPDikLYPFIYHAvESCAETADHSGQRTGT

靶点信息

研究背景	Poly(ADP-ribose) glycohydrolase (PARG) is the major enzyme responsible for the catabolism of poly(ADP-ribose), a reversible covalent-modifier of chromosomal proteins. The protein is found in many tissues and may be subject to proteolysis generating smaller, active products. Several transcript variants encoding different isoforms have been found for this gene.
基因ID	8505

基因名	
Swiss	Q86W56
别名	PARG;PARG99

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn