

ICAM-1/CD54 Rabbit pAb

货号: B10576

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus , Sus scrofa , Rattus norvegicus , Homo sapiens , Anatinae , Homo sapiens Mus musculus IF: Sus scrofa
应用	WB IHC IF/ICC
推荐浓度	WB: 1:1000 - 1:5000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	57kDa
实测分子量	89KDa/92KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,Raji,Mouse lung
细胞定位	Membrane,Single-pass type I membrane protein
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 120-320 of human ICAM-1/CD54 (NP_000192.2).
序列	PLPSWQPVGKNLTLRCQVEGGAPRANLTVLLRGEKELKREPAVGEPAEVTTTVLVRRDHHGANFSCRTELDLRPQGLELF ENTSAPYQLQTFVLPATPPQLVSPRVLEVDVTQGTVCVSLDGLFPVSEAAQVHLALGDQRLNPTVTYGNDSFSAKASVSVTA EDEGTQRLTCAVILGNQSQETLQTVTIYSFPAPNVILTKP

靶点信息

研究背景	This gene encodes a cell surface glycoprotein which is typically expressed on endothelial cells and cells of the immune system. It binds to integrins of type CD11a / CD18, or CD11b / CD18 and is also exploited by Rhinovirus as a receptor.
基因ID	3383
基因名	ICAM1
Swiss	P05362
别名	BB2;CD54;P3.58;ICAM1

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn