

PBEF/Visfatin/NAMPT Rabbit pAb

货号: **B10591**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: human , Mus musculus , Gallus gallus , Rattus norvegicus , Homo sapiens IF: Rattus norvegicus , Homo sapiens IHC: Homo sapiens , Sus scrofa , Mus musculus
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	55kDa
实测分子量	56KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,Mouse liver,Rat liver
细胞定位	Cytoplasm,Nucleus,Secreted
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 181-382 of human PBEF /Visfatin/NAMPT (NP_005737.1).
序列	GNLDGLEYKLHDFGYRGVSSQETAGIGASAHLVNFKGTDTVAGLALIKKYYGTKDPVPGYSVPAAEHSTITAWGKDHEKD AFEHIVTQFSSVPVSVVSDSYDIYNACEKIWGEDLRHLIVSRSTQAPLIIRPDSGNPLDTVLKVLKLEILGKKFPVTENSKGYKLL PPYLRVIQGDGVDINTLQEIVEGMKQKMWSIENIAFGS

靶点信息

研究背景	This gene encodes a protein that catalyzes the condensation of nicotinamide with 5-phosphoribosyl-1-pyrophosphate to yield nicotinamide mononucleotide, one step in the biosynthesis of nicotinamide adenine dinucleotide. The protein belongs to the nicotinic acid phosphoribosyltransferase (NAPRTase) family and is thought to be involved in many important biological processes, including metabolism, stress response and aging. This gene has a pseudogene on chromosome 10.
基因ID	10135
基因名	NAMPT
Swiss	P43490
别名	NAMPT;1110035O14Rik;PBEF;PBEF1;VF;VISFATIN

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn