

BACH1 Rabbit pAb

货号: **B10623**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus , Rattus norvegicus , Homo sapiens IHC: Homo sapiens
应用	WB IP
推荐浓度	WB: 1:500 - 1:1000 IP: 1:500 - 1:1000
理论分子量	81kDa
实测分子量	82KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	K-562
细胞定位	Nucleus
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 410-736 of human BACH1 (NP_001177.1).
序列	TPCQMQLSPAVAKDGEISQKRSECPWLGIRISESPEPGQRTFTTLSSVNC PFISTLSTEGCSSNLEIGNDDYVSEPQQEP CPYACVISLGDDSETDTEGDSESCSAREQECEVKLPFNAQRIISLSRND FQSLLKMHKLTPEQLDCIHDIRRRSKNRIAAQR CRKRKLDICQNLSEIEIKLQSEKESLLKERDHILSTLGETKQNL TGLCQKVCKEAALSQEQIQILAKYSAADCPLSFLISEKDK STPDGELALPSIFSLSDRPPAVLP PCARGNSEPGYARGQESQQMSTATSEQAGPAEQCRQSGGISDFCQMQMTDKCTTDE

靶点信息

研究背景	This gene encodes a transcription factor that belongs to the cap'n'collar type of basic region leucine zipper factor family (CNC-bZip). The encoded protein contains broad complex, tramtrack, bric-a-brac/poxvirus and zinc finger (BTB/POZ) domains, which is atypical of CNC-bZip family members. These BTB/POZ domains facilitate protein-protein interactions and formation of homo- and/or hetero-oligomers. When this encoded protein forms a heterodimer with MafK, it functions as a repressor of Maf recognition element (MARE) and transcription is repressed. Multiple alternatively spliced transcript variants have been identified for this gene.
基因ID	571
基因名	BACH1
Swiss	O14867
别名	BACH1;BACH-1;BTBD24

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn