

CD80 Rabbit pAb

货号: **B10631**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	IF: Mus musculus , Rattus norvegicus FC: Homo sapiens WB: Homo sapiens , Mus musculus
应用	WB IHC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200
理论分子量	18kDa/29kDa/33kDa
实测分子量	50-75KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Rat lung
细胞定位	Membrane,Single-pass type I membrane protein
纯化	Affinity purification

抗原信息

抗原信息	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CD80 (NP_005182.1).
序列	MGHTRRQGTSPSKCPYLNFFQLLVLAGLSHFCSGVIHVTKEVKEVATLSCGHNVSVEELAQTRIWQKEKKMVLTMMSG DMNIWPEYKNRTIFDITNNLS

靶点信息

研究背景	The protein encoded by this gene is a membrane receptor that is activated by the binding of CD28 or CTL A-4. The activated protein induces T-cell proliferation and cytokine production. This protein can act as a receptor for adenovirus subgroup B and may play a role in lupus neuropathy.
------	--

基因ID	941
基因名	CD80
Swiss	P33681
别名	B7;B7-1;B7.1;BB1;CD28LG;CD28LG1;LAB7;CD80

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn