

KIAA1524 Rabbit pAb

货号: **B10634**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens,Mus musculus,Rattus norvegicus Co-IP: Homo sapiens,Mus musculus,Rattus norvegicus
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	84kDa/102kDa
实测分子量	102kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,U-87MG,MCF7,HeLa,A375,Rat testis
细胞定位	Cytoplasm,Membrane,Single-pass membrane protein
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 681-905 of human KIAA 1524 (NP_065941.2).
序列	EVERKNEELSVLLKAQQVESERAQSDIEHLFQHNRKLESVAEEHEILTKSYMELLQRNESTEKKNKDLQITCDSLNKQIETVK KLNESLKEQNEKSIAQLIEKEEQRKEVQNQLVDREHKLANLHQTKVQEEKIKTLQKEREDKEETIDILRKELSRTEQIRKELSI KASSLEVQKAQLEGRLEEKESLVKLQEEELNKHSHMIAMIIHSLSGGKINPETVNLSI

靶点信息

研究背景	Oncoprotein that inhibits PP2A and stabilizes MYC in human malignancies. Promotes anchorage-independ ent cell growth and tumor formation.
基因ID	57650

基因名	KIAA1524
Swiss	Q8TCG1
别名	KIAA1524;CIP2A;p90

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn