

KIFC1 Rabbit pAb

货号: **B10643**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens IHC: Homo sapiens IF: Homo sapiens
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	73kDa
实测分子量	74kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,HT-29,HeLa,A-549
细胞定位	Cytoplasm,Early endosome,Nucleus,centrosome,cytoskeleton,microtubule organizing center,spindle
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human KIFC1 (NP_002254.2).
序列	MDPQRSPLLEVKGNIELKRPLIKAPSQLPLSGSRLKRRPDQMEDGLEPEKKRTRGLGATTKITTSHPRVPSLTTVPQTQGQT TAQKVSKKTGPRCSTAIATGLKNQKVPVAVPVQKSGTSGVPPMAGGKKPSKRPAWDLKGQLCDLNAELKRCRERTQTLD QENQQLQDQLRDAQQVQKALGTERTTLEGLAKVQAQAEQGGQELKNLRACVLELEERLSTQEGLVQELQKKQVELQEE RRGLMSQLEE

靶点信息

研究背景	Minus end-directed microtubule-dependent motor required for bipolar spindle formation. May contribute to movement of early endocytic vesicles (By similarity. Regulates cilium formation and structure (By similarity.
基因ID	3833
基因名	KIFC1
Swiss	Q9BW19
别名	KIFC1;HSET;KNSL2

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn