

GKN1 Rabbit pAb

货号: B10666

产品信息

反应	Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	IHC: Homo sapiens , Mus musculus WB: Mus musculus
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	21kDa
实测分子量	22kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse stomach
细胞定位	Secreted
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 25-185 of human GKN1 (NP_062563.3).
序列	LGVFLAPALANYNINVNDDNNNAGSGQQSVSVNNEHNVANVDNNGWDSWNSIWVDYGNNGFAATRLFQKKTCIVHKM NKEVMPSIQSLDALVKEKKLQGKGPGGPPPKGLMYSVNPKNVDDLSKFGKNIANMCRGIPTYMAEEMQEASLFFYSGTCY TTSVL

靶点信息

研究背景	The protein encoded by this gene is found to be down-regulated in human gastric cancer tissue as compared to normal gastric mucosa.
基因ID	56287

基因名	GKN1
Swiss	Q9NS71
别名	GKN1;AMP18;BRICD1;CA11;FOV;foveolin

产品验证



Western blot analysis of GKN1 expressed in Mouse stomach using GKN1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn