

LAMP1 Rabbit pAb

货号: **B10675**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	IHC: Homo sapiens,Mus musculus WB: Mus musculus
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	38kDa/44kDa
实测分子量	45kDa/120kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse liver
细胞定位	Cell membrane,Endosome membrane,Late endosome,Lysosome membrane,Single-pass type I membrane protein
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 51-270 of human LAMP1 (NP_005552.3).
序列	SVNYDTKSGPKNMTFDLPSDATVVLNRSSCGKENTSDPSLVIAFGRGHTLTlnftrNATRYSVQLMSFVYNLSDTHLFPN ASSKEIKTVESITDIRADIDKKYRCVSGTQVHMNNVTVTlHDAIQAyLSNSSFSRGETRCEQDRPSPTTAPPAPPSPSPSP VPKSPVDKYNVSGTNGTCLlASMGlQLNLTYERKDNTTVTRLLNINPNKTSASGSCG

靶点信息

研究背景	The protein encoded by this gene is a member of a family of membrane glycoproteins. This glycoprotein provides selectins with carbohydrate ligands. It may also play a role in tumor cell metastasis.
------	---

基因ID	3916
基因名	LAMP1
Swiss	P11279
别名	LAMP1;CD107a;LAMPa;LGP120

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn