

KIF5A Rabbit pAb

货号: B10699

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	IF: Mus musculus
应用	WB IF/ICC
推荐浓度	WB: 1:1000 - 1:3000 IF/ICC: 1:50 - 1:200
理论分子量	117kDa
实测分子量	129kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Rat brain
细胞定位	Cytoplasm,cytoskeleton,perinuclear region
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 933-1032 of human KIF5A (NP_004975.2).
序列	TNPYGTRSPECISYTNLSLFQNYQNLYLQATPSSTSDMYFANSCTSSGATSSGGPLASYQKANMDNGNATDINDNRSDDLPCGYEAEDQAKLFPLHQETAAS

靶点信息

研究背景	This gene encodes a member of the kinesin family of proteins. Members of this family are part of a multisubunit complex that functions as a microtubule motor in intracellular organelle transport. Mutations in this gene cause autosomal dominant spastic paraplegia 10.
基因ID	3798

基因名	KIF5A
Swiss	Q12840
别名	KIF5A;D12S1889;MY050;NEIMY;NKHC;SPG10

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn