

CDK19 Rabbit pAb

货号: B10713

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus
应用	WB IHC
推荐浓度	WB: 1:100 - 1:500 IHC: 1:50 - 1:200
理论分子量	49kDa/56kDa
实测分子量	57KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	DU145,Mouse brain
细胞定位	cytosol,nucleoplasm,nucleus,perinuclear region of cytoplasm
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 250-350 of human CDK19 (XP_005266928.1).
序列	QQQNQHQQPTAPPQQAAPPQAPPPQQNSTQTNGTAGGAGAGVGGTGAGLQHSQDSSLNQVPPNKKPRLGPSGANS GGPVMPSDQYHSSRLNYQSSVQGS

靶点信息

研究背景	This gene encodes a protein that is one of the components of the Mediator co-activator complex. The Mediator complex is a multi-protein complex required for transcriptional activation by DNA binding transcription factors of genes transcribed by RNA polymerase II. The protein encoded by this gene is similar to cyclin-dependent kinase 8 which can also be a component of the Mediator complex. Alternative splicing results in multiple transcript variants encoding distinct isoforms.
------	--

基因ID	23097
基因名	CDK19
Swiss	Q9BWU1
别名	CDK19;CDC2L6;CDK11;bA346C16.3

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn