

KTN1 Rabbit pAb

货号: B10730

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Arabidopsis thaliana
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	149kDa/150kDa/152kDa/156kDa
实测分子量	178kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	K-562,SKOV3,Raji,HeLa,HT-1080,Mouse heart,Rat testis
细胞定位	Endoplasmic reticulum membrane,Single-pass type II membrane protein
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human KTN1 (NP_001072989.1).
序列	MEFYESAYFIVLIPSIVITVIFLFFWLFMKETLYDEVLAQKQREQKLIPTKTDKKKAEEKKNKKKEIQNGNLHESDSSESVPRDFKLSDALAVEDDQVAPV

靶点信息

研究背景	This gene encodes an integral membrane protein that is a member of the kinectin protein family. The encoded protein is primarily localized to the endoplasmic reticulum membrane. This protein binds kinesin and may be involved in intracellular organelle motility. This protein also binds translation elongation factor-delta and may be involved in the assembly of the elongation factor-1 complex. Alternate splicing results in multiple transcript variants of this gene.
------	--

基因ID	3895
基因名	KTN1
Swiss	Q86UP2
别名	KTN1;CG1;KNT;MU-RMS-40.19;kinectin

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn