

PPT2 Rabbit pAb

货号: B10732

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	IHC: Homo sapiens
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	34kDa
实测分子量	36kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse liver
细胞定位	Lysosome
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 28-302 of human PPT2 (NP_005146.4).
序列	APAPHRASYKPVIVVHGLFDSSYSFRHLLEYINETHPGTVVTVLDLFDGRESLRPLWEQVQGFREAVVPIMAKAPQGVHLIC YSQGGLVCRALLSVMDHNVDSFISLSSPQMGOYGDYDLKWLFPSTMRSNLYRICYSPWQGEFSICNYWHDPHHDDL YLNASSFLALINGERDHPNATVWRKNFLRVGHLVLIGGPDDGVITPWQSSFFGFYDANETVLEMEEQLVYLDRDSFGLKTL LARGAIVRCPMAGISHTAWHSNRTLYETCIEPWLS

靶点信息

研究背景	This gene encodes a member of the palmitoyl-protein thioesterase family. The encoded glycosylated lyso somal protein has palmitoyl-CoA hydrolase activity in vitro, but does not hydrolyze palmitate from cystein e residues in proteins. Alternative splicing results in multiple transcript variants. Read-through transcripti on also exists between this gene and the downstream EGFL8 (EGF-like-domain, multiple 8) gene.
------	--

基因ID	9374
基因名	PPT2
Swiss	Q9UMR5
别名	PPT2;C6orf8;G14;PPT-2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn