

SLC14A1 Rabbit pAb

货号: B10746

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:100 - 1:500 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	42kDa/48kDa
实测分子量	35-45KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	RT4
细胞定位	basolateral plasma membrane,plasma membrane
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 100 to the C-terminus of human SLC 14A1 (NP_001295208.1).
序列	QIYGCDNPWTGGIFLGAILLSSPLMCLHAAIGSLLGIAAGLSLSAPFEDIYFGLWGFNSSLACIAMGGMFALTWQTHLLALGCALFTAYLGVGMANFMAEVLGPACTWPFCLATLLFLIMTTKNSNIYKMPLSKVTYPEENRIFYLQAKKRMVESPL

靶点信息

研究背景	The protein encoded by this gene is a membrane transporter that mediates urea transport in erythrocyte s. This gene forms the basis for the Kidd blood group system.
基因ID	6563

基因名	SLC14A1
Swiss	Q13336
别名	SLC14A1;HUT11;HsT1341;JK;RACH1;RACH2;UT-B1;UT1;UTE

产品验证



Western blot analysis of SLC14A1 expressed in RT4 using SLC14A1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn