

NOL3 Rabbit pAb

货号: B10779

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:10 - 1:100
理论分子量	22kDa/24kDa/29kDa
实测分子量	26kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	BT-474,MCF7,SW620,U-937,SW480,Mouse heart,Mouse lung,Mouse skeletal muscle,Rat heart
细胞定位	Cytoplasm,Lipid-anchor,Membrane,Mitochondrion,Nucleus,Sarcoplasmic reticulum,nucleolus
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-208 of human NOL3 (NP_003937.1).
序列	MGNAQERPSETIDRERKRLVETLQADSGLLLDALLARGVLTGPEYEALDALPDAERRVRRLLLLVQKGGEAACQELLRCAQ RTAGAPDPAWDWQHVGPGYRDRSYDPPCPGHWTPPEAPGSGTTCPGLPRASDPDEAGGPEGSEAVQSGTPPEEPEPELE AEASKEAPEPEPEPELEPEAEAEPEPELEPEPDPEPEPDFEERDESEDS

靶点信息

研究背景	This gene encodes an anti-apoptotic protein that has been shown to down-regulate the enzyme activities of caspase 2, caspase 8 and tumor protein p53. Multiple transcript variants encoding different isoforms have been found for this gene.
------	---

基因ID	8996
基因名	NOL3
Swiss	O60936
别名	NOL3;ARC;FCM;MYP;NOP;NOP30

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn