

YIPF7 Rabbit pAb

货号: B10799

产品信息

反应	Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	7kDa/14kDa/23kDa/30kDa
实测分子量	31kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Rat heart,Rat liver
细胞定位	Endoplasmic reticulum membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-80 of human YIPF7 (NP_872398.2).
序列	MDLLKISHTKLHLLLEDLSIKNKQRMSNLAQFDSDFYQSNFTIDNQEQSGNDSNAYGNLYGSRKQQAGEQPQPASFVPSE M

靶点信息

研究背景	Predicted to be involved in endoplasmic reticulum to Golgi vesicle-mediated transport and vesicle fusion with Golgi apparatus. Predicted to be located in endoplasmic reticulum membrane. Predicted to be integral component of membrane. Predicted to be active in trans-Golgi network.
基因ID	285525

基因名	YIPF7
Swiss	Q8N8F6
别名	YIPF7;FinGER9

产品验证



Western blot analysis of YIPF7 expressed in Rat heart,Rat liver using YIPF7 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn