

SAMD3 Rabbit pAb

货号: B10803

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	25kDa/61kDa/63kDa
实测分子量	61kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,Mouse brain,Rat brain
细胞定位	
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human SAMD3 (NP_001017373.2).
序列	METWSVEQVCSWLVEKNLGELVHRFQEEEVSGAALLALNDRMVQQLVKKIGHQAVLMDLIKKYQNTQGLKSPENPKKA ALVMQTEAARDYRDEESSSPARHGEQMPSFYPAENLDNGLIDQRVLKQRRNVKQILARSKALQWTKSYVLPEFPYDVKC MLAEQKCPDHSMRIRIIEFLQADMTKYLEGSLYPSTQQYNDVVNALLQAHPF

靶点信息

研究背景	SAMD3 (Sterile Alpha Motif Domain Containing 3) is a Protein Coding gene. Diseases associated with SAM D3 include Diaphyseal Medullary Stenosis With Malignant Fibrous Histiocytoma. An important paralog of t his gene is CNKSR2.
基因ID	154075

基因名	SAMD3
Swiss	Q8N6K7
别名	SAMD3

产品验证



Western blot analysis of SAMD3 expressed in U-251MG,Mouse brain,Rat brain using SAMD3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn