

ZNF644 Rabbit pAb

货号: B10817

产品信息

反应	Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	11kDa/139kDa/149kDa
实测分子量	150kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse liver
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human ZNF644 (NP_958357.1).
序列	MRSFLQQDVNKTKSRLNVLNGLANNMDDLKINTDITGAKEELLDDNNFISDKESGVHKPKDCQTSFQKNNTLTLPEELSK DKSENALSGGQSSLFIHAGAPTVSSENFILPKGAAVNGPVSHSSLTKTSNMNKGSVSLTTGQPVDQPTTESCSTLKVAAD LQLSTPQKASQHQLFLFLSDVAHAKNPHTSNKKLPTSASVGCDIQNSVGSNIKSDGTLINQVEVGEDGEDLLVKDDCVNT VTGISSGTDGFRSENDTNWDPQKEFIQLMTNEETVDKAPPHSKIGLEKKRKRKMDVSKI

靶点信息

研究背景	The protein encoded by this gene is a zinc finger transcription factor that may play a role in eye develop ment. Defects in this gene have been associated with high myopia. Three transcript variants encoding tw o different isoforms have been found for this gene.
基因ID	84146

基因名	ZNF644
Swiss	Q9H582
别名	ZNF644;BM-005;MYP21;NatF;ZEP-2

产品验证



Western blot analysis of ZNF644 expressed in Mouse liver using ZNF644 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn