

# NLGN1 Rabbit pAb

货号: B10848

## 产品信息

|       |                                                                                                                         |
|-------|-------------------------------------------------------------------------------------------------------------------------|
| 反应    | Mouse                                                                                                                   |
| 宿主    | Rabbit                                                                                                                  |
| 克隆性   | Polyclonal                                                                                                              |
| 预测反应  |                                                                                                                         |
| 应用    | WB                                                                                                                      |
| 推荐浓度  | WB: 1:500 - 1:2000                                                                                                      |
| 理论分子量 | 92kDa/93kDa                                                                                                             |
| 实测分子量 | 105kDa                                                                                                                  |
| 形式    | Liquid                                                                                                                  |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.                    |
| 偶联物   | Unconjugated                                                                                                            |
| 阳性对照  | Mouse kidney                                                                                                            |
| 细胞定位  | Cell junction,Cell membrane,Single-pass type I membrane protein,postsynaptic cell membrane,postsynaptic density,synapse |
| 纯化    | Affinity purification                                                                                                   |

## 抗原信息


|      |                                                                                                          |
|------|----------------------------------------------------------------------------------------------------------|
| 抗原信息 | A synthetic peptide corresponding to a sequence within amino acids 200-300 of human NLGN1 (NP_055747.1). |
| 序列   | YGNVIVITVNYRLGVLGFLSTGDQAAKGNYGLLDLIQALRWTSENIGFFGGDPLRITVFGSGAGGSCVNLLTSHYSEGNRWSNSTKGLFQRAIAQSGTAL     |

## 靶点信息

|      |                                                                                                                                                                                                                                                     |
|------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This gene encodes a member of a family of neuronal cell surface proteins. Members of this family may act as splice site-specific ligands for beta-neurexins and may be involved in the formation and remodeling of central nervous system synapses. |
| 基因ID | 22871                                                                                                                                                                                                                                               |

|       |           |
|-------|-----------|
| 基因名   | NLGN1     |
| Swiss | Q8N2Q7    |
| 别名    | NLGN1;NL1 |

产品验证



Western blot analysis of NLGN1 expressed in Mouse kidney using NLGN1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)