

SYCP2 Rabbit pAb

货号: B10850

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	175kDa
实测分子量	176kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,Mouse testis,Rat testis
细胞定位	Chromosome,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1231-1530 of human SYCP2 (NP_055073.2).
序列	EIESPHINENYIQSKREESHCLASSLSKSSEGREKTWFDMPCDATHVSGPTQHLSRKRIYIEDNLSNSNEVEMEKEGERRANL LPKKLKIEDADHHIHKMSESVSLSSTNDFSIPWETWQNEFAGIEMTYETYERLNSEFKRRNNIRHKMLSFTTQSWKTAQ QHLRTMNHQSQDSRIKKLDKFQFIIIELENFEKDSQSLKDLEKEFVDFWEKIFQKFSAYQKSEQQLHLLKTSLAKSVFCN TDSEETVFTSEMCLMKEDMKVLQDRLLKDMLEEELLNVRRELMVSFMSHERNANV

靶点信息

研究背景	The synaptonemal complex is a proteinaceous structure that links homologous chromosomes during the prophase of meiosis. The protein encoded by this gene is a major component of the synaptonemal complex and may bind DNA at scaffold attachment regions. The encoded protein requires synaptonemal complex protein 3, but not 1, for inclusion in the synaptonemal complex.
------	---

基因ID	10388
基因名	SYCP2
Swiss	Q9BX26
别名	SYCP2;SCP-2;SCP2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn