

GPRASP1 Rabbit pAb

货号: B10854

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	156kDa
实测分子量	157kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Rat eye,Rat testis
细胞定位	Cytoplasm
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 200-500 of human GPR ASP1 (NP_055525.3).
序列	SSLFWSGDEVTAKFHPGNRVKDSNRSMHMANQEANTMSRSQTNQELYIASSSGSEDESVKTPWFFWARDKTNTWSGP REDPNRSRFRSKKEVYVESSSGSEHEDHLESWFGAGKEAKFRSKMRAGKEANNRARHRAKREACIDFMPGSIDVIKKES CFWPEENANTFSRPMIKKEARARAMTKEEAKTKARARAKQEARSEEEALIGTWFWATDESSMADEASIESSLQVEDESIIG SWFWTEEEASMG TGASSKSRPRTDGERIGDSLFGAREKTSMTKGAEATSESILAADDEQVIIGS

靶点信息

研究背景	This gene encodes a member of the GPRASP (G protein-coupled receptor associated sorting protein) family. The protein may modulate lysosomal sorting and functional down-regulation of a variety of G-protein coupled receptors. It targets receptors for degradation in lysosomes. The receptors interacting with this sorting protein include D2 dopamine receptor (DRD2), delta opioid receptor (OPRD1), beta-2 adrenergic receptor (ADRB2), D4 dopamine receptor (DRD4) and cannabinoid 1 receptor (CB1R). Multiple alternatively spliced transcript variants encoding the same protein have been identified.
基因ID	9737
基因名	GPRASP1
Swiss	Q5JY77
别名	GPRASP1;GASP;GASP-1;GASP1

产品验证



Western blot analysis of GPRASP1 expressed in Rat eye,Rat testis using GPRASP1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn