

KCNK6 Rabbit pAb

货号: B10858

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	19kDa/33kDa
实测分子量	36kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	BxPC-3,HT-29,K-562,A-549,293T,Mouse pancreas,Mouse heart
细胞定位	Membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 254-313 of human KCN K6 (NP_004814.1).
序列	QTFRHHVSDLHGLTEILLPPPCPASFNADEDDRVDILGPQESHQQLSASSHTDYASIPR

靶点信息

研究背景	This gene encodes one of the members of the superfamily of potassium channel proteins containing two pore-forming P domains. This channel protein, considered an open rectifier, is widely expressed. It is stimulated by arachidonic acid, and inhibited by internal acidification and volatile anaesthetics.
基因ID	9424
基因名	KCNK6

Swiss	Q9Y257
别名	KCNK6;K2p6.1;KCNK8;TOSS;TWIK-2;TWIK2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn