

MKRN3 Rabbit pAb

货号: B10864

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	55kDa
实测分子量	56kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse kidney
细胞定位	
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human MKRN3 (NP_005655.1).
序列	MEEPAAPSEAHEAAGAQAAGAEAAAREGVSGPDLPVCEPSGESAAPDSALPHAARGWAPFPVAPVPAHLRRGGGLRPAPASG GGAWPSPLPSRSSGIWTKQIICRYIYHGQCKEGENCRYSHDLSGRKMATGGVSPPGASAGGGPSTAAHIEPPTQEVAEA PPAASSLSLPVIGSAAERGFFEAERDNADRGAAAGGAGVESWADAIEFVPGQPYPYRGRWVASAPEPLQSSETERKQMAVG SGLRFCYYASRGVCFRGESCMY

靶点信息

研究背景	The protein encoded by this gene contains a RING (C3HC4) zinc finger motif and several C3H zinc finger motifs. This gene is intronless and imprinted, with expression only from the paternal allele. Disruption of the imprinting at this locus may contribute to Prader-Willi syndrome. An antisense RNA of unknown function has been found overlapping this gene.
------	---

基因ID	7681
基因名	MKRN3
Swiss	Q13064
别名	MKRN3;CPPB2;D15S9;RNF63;ZFP127;ZNF127

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn