

PLXNB3 Rabbit pAb

货号: B10865

产品信息

反应	Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	206kDa/209kDa
实测分子量	170kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain
细胞定位	Cell membrane,Single-pass type I membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1050-1250 of human PLXNB3 (NP_001156729.1).
序列	RGGGRLIRVRGTGLDVVQRPLLSVWLEADAEVQASRAQPQDPQPRRSCGAPAADPQACIQLGGGLLQCSTVCSVNSSSL LLCRSPAVPDRAHPQRVFFTLDNVQVDFASASGGQGFLYQPNRLAPLSREGPARPYRLKPGHVLDVEGEGLNLGISKEE VRVHIGRGECLVKTLTRTHLYCEPPAHAPQPANGSGLPQFVV

靶点信息

研究背景	The protein encoded by this gene is a member of the plexin family. It functions as a receptor for semaphorin 5A, and plays a role in axon guidance, invasive growth and cell migration. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
基因ID	5365

基因名	PLXNB3
Swiss	Q9ULL4
别名	PLXNB3;PLEXB3;PLEXR;PLXN6;plexin-B3

产品验证



Western blot analysis of PLXNB3 expressed in Mouse brain using PLXNB3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn