

FBLN1 Rabbit pAb

货号: **B10872**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	61kDa/65kDa/74kDa/77kDa
实测分子量	77KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,A-549,U-87MG,Mouse brain,Rat brain,Rat heart
细胞定位	Secreted,extracellular matrix,extracellular space
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 570-670 of human FBLN 1 (NP_001987.2).
序列	RLPCHENRECSKLPLRITYYHLSFPTNIQAPAVVFRMGPPSSAVPGDSMQLAITGGNEEGFFTRKVSPhSGVVAlTKPVPEPRDLLLLTVKMDLSRHGTvSS

靶点信息

研究背景	Fibulin 1 is a secreted glycoprotein that becomes incorporated into a fibrillar extracellular matrix. Calciu m-binding is apparently required to mediate its binding to laminin and nidogen. It mediates platelet adhe sion via binding fibrinogen. Four splice variants which differ in the 3' end have been identified. Each varia nt encodes a different isoform, but no functional distinctions have been identified among the four variant s.
基因ID	2192

基因名	FBLN1
Swiss	P23142
别名	FBLN1;FBLN;FIBL1;fibulin-1

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn