

ETFB Rabbit pAb

货号: B10878

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:100
理论分子量	27kDa/37kDa
实测分子量	30kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,A-549,293T,BxPC-3,Mouse brain,Mouse heart,Mouse kidney,Mouse liver
细胞定位	Mitochondrion matrix
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 80-230 of human ETFB (NP_001976.1).
序列	AMGADRGIHVEVPPAEERLGPLQVARVLAKLAEKEKVDLVLLGKQAIDDDCNQTGQMTAGFLDWPQGTFASQVTLEG DKLKVEREIDGGLETLRCLKPAVVTTADLRRLNEPRYATLPNIMKAKKKKIEVIKPGDLGVDLTSKLSVISVEDP

靶点信息

研究背景	This gene encodes electron-transfer-flavoprotein, beta polypeptide, which shuttles electrons between primary flavoprotein dehydrogenases involved in mitochondrial fatty acid and amino acid catabolism and the membrane-bound electron transfer flavoprotein ubiquinone oxidoreductase. The gene deficiencies have been implicated in type II glutaricaciduria. Alternatively spliced transcript variants have been found for this gene.
------	---

基因ID	2109
基因名	ETFB
Swiss	P38117
别名	ETFB;FP585;MADD

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn