

RNF10 Rabbit pAb

货号: B10907

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:100 - 1:500
理论分子量	89kDa/90kDa
实测分子量	90KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Rat testis
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 552-811 of human RNF10 (NP_055683.3).
序列	SMSQDVQRHRYLSHLPLTCEFSICELALQPPVSKETLEMFSDDIEKRKRQRQKKAREERRRERRIEIEENKKQGYEVHILENLQQFPFNSYTCSSDSALGPTSTEGHGALSIPLSRSPGSHADFLLTPLSPTASQGSPSFCVGSLEEDSPFPSFAQMLRVGKAKADVWPKTAPKKDENSELVPPAPVDSGSDNSDRVPVPSFQNSFSQAIEAAFMKLDTPATSDPLSEEKGGKKRKQKQKLLFSTSVVHTK

靶点信息

研究背景	The protein encoded by this gene contains a ring finger motif, which is known to be involved in protein-protein interactions. The specific function of this protein has not yet been determined. EST data suggests the existence of multiple alternatively spliced transcript variants, however, their full length nature is not known.
------	---

基因ID	9921
基因名	RNF10
Swiss	Q8N5U6
别名	RNF10;RIE2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn