

PLA2G4C Rabbit pAb

货号: B10954

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	IF/ICC
推荐浓度	IF/ICC: 1:50 - 1:200
理论分子量	59kDa/60kDa/62kDa
实测分子量	59kDa/60kDa/62kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Lipid-anchor,Membrane
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-189 of human PLA2G4 C (NP_003697.2).
序列	MGSSEVSIIPGLQKEEKAAVERRRRLHVLKALKKLRIEADEAPVVAVLGSGGGLRAHIACLGVLSEMKEQGLLDAVTYLAGVSGSTWAISSLYTNDGDMEALEADLKHFRFTRQEWDLAKSLQKTIQAARSENYSLTDFWAYMVISKQTRELPESHLSNMKKPV EEGTLPYPIFAAIDNDLQPSWQEARAP

靶点信息

研究背景	This gene encodes a protein which is a member of the phospholipase A2 enzyme family which hydrolyzes glycerophospholipids to produce free fatty acids and lysophospholipids, both of which serve as precursor s in the production of signaling molecules. The encoded protein has been shown to be a calcium-indepen dent and membrane bound enzyme. Multiple transcript variants encoding different isoforms have been fo und for this gene.
------	---

基因ID	8605
基因名	PLA2G4C
Swiss	Q9UP65
别名	PLA2G4C;CPLA2-gamma

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn