

SOCS3 Rabbit pAb

货号: B10962

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: PK15 , Homo sapiens , Mus musculus , Rattus norvegicus IHC: Mouse liver , Rattus norvegicus
应用	WB IHC IF/ICC
推荐浓度	WB: 1:100 - 1:500 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	24kDa
实测分子量	30KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	RAW264.7
细胞定位	cytosol
纯化	Affinity purification


抗原信息

抗原信息	A synthetic peptide corresponding to a sequence within amino acids 100 to the C-terminus of human SOC S3 (NP_003946.3).
序列	GSFSLQSDPRSTQPVPRFDCVLKLVHHYMPPPGAPSFPSPTPEPSSEVPEQPSAQPLPGSPPRRAYYYISGGEKIPLVLSRP LSSNVATLQHLCRKTVNGHLDSYEKVTQLPGPIREFLDQYDAPL

靶点信息

研究背景	This gene encodes a member of the STAT-induced STAT inhibitor (SSI), also known as suppressor of cytokine signaling (SOCS), family. SSI family members are cytokine-inducible negative regulators of cytokine signaling. The expression of this gene is induced by various cytokines, including IL6, IL10, and interferon (IFN)-gamma. The protein encoded by this gene can bind to JAK2 kinase, and inhibit the activity of JAK2 kinase. Studies of the mouse counterpart of this gene suggested the roles of this gene in the negative regulation of fetal liver hematopoiesis, and placental development.
基因ID	9021
基因名	SOCS3
Swiss	O14543
别名	SOCS3;ATOD4;CIS3;Cish3;SOCS-3;SSI-3;SSI3

产品验证



Western blot analysis of SOCS3 expressed in RAW264.7 using SOCS3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn