

IKKε Rabbit pAb

货号: B10965

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: 293T,Flag-IKKε,THP-1,Hela,BMDM , Homo sapiens , Mus musculus , Sus scrofa IP: Homo sapiens
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	70kDa/80kDa
实测分子量	80KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SW480,DU145,MCF7,Mouse liver
细胞定位	Cytoplasm,Nucleus,PML body
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 497-716 of human IKKε (NP_054721.1).
序列	SRLRTLAEVLSRCSQNITETQESLSSLNRELVKSRDQVHEDRSIQIQCCLDKMNFIYKQFKKSRMRPGLGYNEEQIHKLDK VNFSLAKRLLQVQEECVQKYQASLVTHGKRMRVVHETRNHLRLVGCSVAACNTEAQGVQESLSKLEELSHQLQDRA KGAQASPPPIAPYPSPTRKDLLLLHMQELCEGMKLLASDLLDNNRIERLNRVPAPPDV

靶点信息

研究背景	IKBKE is a noncanonical I-kappa-B (see MIM 164008) kinase (IKK) that is essential for regulating antiviral signaling pathways. IKBKE has also been identified as a breast cancer (MIM 114480) oncogene and is amplified and overexpressed in over 30% of breast carcinomas and breast cancer cell lines (Hutti et al., 2009 [PubMed 19481526]).[supplied by OMIM, Oct 2009]
基因ID	9641
基因名	IKBKE
Swiss	Q14164
别名	IKBKE;IKK-E;IKK-i;IKKE;IKKI

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn