

SOD1 Rabbit pAb

货号: B11081

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: 293T , Mus musculus , H. cordata , Homo sapiens , Rattus norvegicus , Bos taurus , Eriocheir sinensis H.Milne-Edwards Co-IP: Mus musculus IHC: Rattus norvegicus
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	15kDa
实测分子量	16-18KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,293T,MCF7,Mouse liver,Mouse brain,Rat lung,Rat liver,Rat brain
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 2-154 of human SOD1 (NP_000445.1).
序列	ATKAVCVLKGDGPVQGIINFEQKESNGPVKVWGSIKGLTEGLHGFHVHEFGDNTAGCTSAGPHFNPLSRKHGGPKDEER HVGD LGNV TADK DG VADVSIEDSVISL SGDHCIIGRTL VVHEKADDLGKGGNEESTKTGNAGSRLACGVIGIAQ

靶点信息

研究背景	The protein encoded by this gene binds copper and zinc ions and is one of two isozymes responsible for destroying free superoxide radicals in the body. The encoded isozyme is a soluble cytoplasmic protein, acting as a homodimer to convert naturally-occurring but harmful superoxide radicals to molecular oxygen and hydrogen peroxide. The other isozyme is a mitochondrial protein. Mutations in this gene have been implicated as causes of familial amyotrophic lateral sclerosis. Rare transcript variants have been reported for this gene.
基因ID	6647
基因名	SOD1
Swiss	P00441
别名	SOD1;ALS;ALS1;HEL-S-44;IPOA;SOD;hSod1;homodimer

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn