

GSDMD (Full Length+N terminal) Rabbit pAb

货号: B11083

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens , Rattus norvegicus , Mus musculus , fish IHC: Mus musculus IF: Mus musculus , Bos taurus , Rattus norvegicus WB、IHC: Mus musculus
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	53kDa
实测分子量	53kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	cytosol,extracellular region,extracellular space,NLRP3 inflammasome complex,nucleoplasm,plasma membrane
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-275 of human GSDMD (NP_079012.3).
序列	MGSAFERVVRVQELDHGGEFIPVTSLSSTGFQPYCLVVRKPSSSWFWKPRYKCVNLSIKDILEPDAAEPDVQRGRSF HFYDAMDGQIQGSVELAAPGQAKIAGGAAVSDSSSTSMNVYSLSDPNTWQTLLHERHLRQPEHKVLQQLRSRGDNVY VVTEVLQTQKEVEVTRTHKREGSGRFSLPGATCLQGEGQGHLSQKKTVTIPSGSTLAFRVAQLVIDSDLVDVLLFPDKKQRT FQPPATGHKRSTSEGAWPQLPSGLSMMRCLHNFLTD

靶点信息

研究背景	Gasdermin D is a member of the gasdermin family. Members of this family appear to play a role in regulation of epithelial proliferation. Gasdermin D has been suggested to act as a tumor suppressor. Alternatively spliced transcript variants have been described.
基因ID	79792
基因名	GSDMD
Swiss	P57764
别名	DF5L;DFNA5L;FKSG10;GSDMDC1;GSDMD

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn