

Slug Rabbit pAb

货号: B11101

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: mouse cells , MDA-MB-231 , Mouse kidney cortex , Mus musculus , Homo sapiens , Homo sapiens、Mus musculus IHC: Homo sapiens IF: Homo sapiens
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	29kDa
实测分子量	30KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Mouse spleen
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-80 of human Slug (NP_003059.1).
序列	MPRSFLVKKHFNASKKPNYSELDTHTVIISPLYESYSMPVIPQPEILSSGAYSPITVWTTAAPFHAQLPNGLSPLSGYS

靶点信息

研究背景	This gene encodes a member of the Snail family of C2H2-type zinc finger transcription factors. The encoded protein acts as a transcriptional repressor that binds to E-box motifs and is also likely to repress E-cadherin transcription in breast carcinoma. This protein is involved in epithelial-mesenchymal transitions and has antiapoptotic activity. Mutations in this gene may be associated with sporadic cases of neural tube defects.
------	---

基因ID	6591
基因名	SNAI2
Swiss	O43623
别名	SNAI2;SLUG;SLUGH;SLUGH1;SNAIL2;WS2D

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn