

# Survivin Rabbit pAb

货号: B11126

产品信息

|       |                                                                                                                                                                                     |
|-------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                                                                                                                     |
| 宿主    | Rabbit                                                                                                                                                                              |
| 克隆性   | Polyclonal                                                                                                                                                                          |
| 预测反应  | <b>WB:</b> Mus musculus , Homo sapiens , Rattus norvegicus , African green monkey<br><b>IF:</b> Mus musculus , Homo sapiens<br><b>ChIP:</b> Mus musculus<br><b>FC:</b> Homo sapiens |
| 应用    | WB IHC                                                                                                                                                                              |
| 推荐浓度  | <b>WB:</b> 1:100 - 1:500<br><b>IHC:</b> 1:50 - 1:200                                                                                                                                |
| 理论分子量 | 8kDa/13kDa/15kDa/16kDa/18kDa                                                                                                                                                        |
| 实测分子量 | 16KDa                                                                                                                                                                               |
| 形式    | Liquid                                                                                                                                                                              |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.                                                                                |
| 偶联物   | Unconjugated                                                                                                                                                                        |
| 阳性对照  | HeLa,293F,Jurkat                                                                                                                                                                    |
| 细胞定位  | Chromosome,Cytoplasm,Midbody,Nucleus,centromere,cytoskeleton,kinetochore,spindle                                                                                                    |
| 纯化    | Affinity purification                                                                                                                                                               |


抗原信息

|      |                                                                                                                                                 |
|------|-------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 1-142 of human Survivin (NP_001159.2).                            |
| 序列   | MGAPTLPPAWQPFLKDHRISTFKNWPFLEGCACTPERMAEAGFIHCPTENEPDLAQCFKCFKELEGWEPDDDDPIEEHKKHSSGCAFLSVKKQFEELTLGEFLKLDREKAKNIAKETNNKKKEFEETAEKVRRRAIEQLAAMD |

靶点信息

|       |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
|-------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景  | This gene is a member of the inhibitor of apoptosis (IAP) gene family, which encode negative regulatory proteins that prevent apoptotic cell death. IAP family members usually contain multiple baculovirus IAP repeat (BIR) domains, but this gene encodes proteins with only a single BIR domain. The encoded proteins also lack a C-terminus RING finger domain. Gene expression is high during fetal development and in most tumors, yet low in adult tissues. Alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. |
| 基因ID  | 332                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| 基因名   | BIRC5                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Swiss | O15392                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| 别名    | API4;EPR-1;BIRC5;Survivin                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)